



## XML Export Manual:


1. Open the Dysmorphology Access DB Input Tool (the switch board should load automatically)




2. Locate and click on the “Export XML files” button




3. A “save-as” dialog box should be brought up




4. Name the file anything you want (in this example I am name it “LRB-data”)



5. Notice I am placing it on the desktop (make sure you are able to find the files created)



6. Four files are created for you (LRB-data.htm, LRB-data.xml, LRB-data.xsd, and LRB-data.xsl)



7. The only file we want upload is the **LRB-data.xml**

